

TECHNICAL DATA SHEET



QUBIX™ PDX

Monofilament Polydioxanone Suture

This versatile suture is supplied in pre-cut lengths, available in USP sizes 4-0 to 1 (multisized), with needles of various types and shapes, and needle lengths ranging from 20 to 50 mm.

Technical Specifications

Product Codes	Description
QN 9210	Size 0, 70 cm, 1/2 Circle Round Body 30 mm
QN 9352	Size 1, 90 cm, 1/2 Circle Round Body Heavy 40 mm
QN 9826	Size 1, 70 cm, 1/2 Circle Reverse Cutting 23 mm
QN 9221	Size 1, 90 cm, 1/2 Circle Reverse Cutting OR8 Heavy 40 mm
QN 9261	Size 0 loop, 150 cm, 1/2 Circle Round Body Heavy 50 mm
QN 9262	Size 1 loop, 150 cm, 1/2 Circle Round Body Heavy 50 mm
QN 9262S	Size 1 loop, 110 cm, 1/2 Circle Round Body Heavy 50 mm
QN 9133	Size 2-0, 70 cm, 1/2 Circle Round Body 30 mm
QN 9132	Size 3-0, 70 cm, 1/2 Circle Round Body 30 mm
QN 9237	Size 3-0, 70 cm, 1/2 Circle Round Body 20 mm
QN 9304	Size 4-0, 70 cm, 1/2 Circle Round Body 20 mm

Sterility: Ethylene Oxide

Packaging: Box of 12 units

Usage: Single-use only