

## TECHNICAL DATA SHEET



QUBIX™ POLYESTER

*Braided and Coated Polyester Suture*

The suture is supplied in pre-cut lengths, available in USP sizes 4-0 to 5, either non-neededled or with needle to various needle types and shapes, with needle lengths ranging from 16–50 mm.

### Technical Specifications

Product Codes	Description
QN 672	Size 4-0, 76cm 3/8circle RB 16mm
QN 632	Size 3-0, 76cm 1/2 circle RB 25mm
QN 617	Size 2-0, 90cm 1/2 C TCDA 17mm
QN 677	Size 2-0, 76cm 1/2 C TCDA 25mm
QN 679	Size 0, 100cm 1/2 Circle RC 50mm
QN 680	Size 1, 100cm 1/2 Circle RC 50mm
QN 643	Size 2, 100cm 1/2 Circle TCH 45mm
QN 645	Size 5, 100cm 1/2 Circle Trocar point Heavy 55mm
QN 646	Size 5, 75cm 1/2 Circle cutting 50mm
QN 4842	Size 2, 75cm 1/2 Circle Taper End Cutting Heavy 45mm
QN 4843	Size 5, 75cm 1/2 Circle Taper End Cutting Heavy 45mm
QN 4845	Size 5, 75cm 1/2 Circle Taper End Cutting Double Armed 55mm

**Sterility:** Ethylene Oxide

**Packaging:** Superior PP box packaging (12 units)

**Usage:** Single-use only