

TECHNICAL DATA SHEET

QUBIX™ SILK

Black Braided Silk Suture

The suture is supplied in pre-cut lengths, available in USP sizes 5-0 to 4, either non-needled or attached to various needle types and shapes, with needle lengths ranging from 12–60 mm.

Technical Specifications

Product Codes	Description
QN 5080	Size 5-0, 76cm, 3/8 Circle RC 12 mm
QN 5082	Size 4-0, 76cm, 3/8 circle RB 16mm
QN 5049	Size 4-0, 76cm, 1/2 circle RB 16mm
QN 5001	Size 4-0, 76cm, 3/8 circle cutting 16mm
QN 5086	Size 4-0, 76cm, 1/2 circle RB 20mm
QN 5002	Size 3-0, 76cm, 3/8 circle cutting 16mm
QN 5087	Size 3-0, 76cm, 1/2 circle RB 20mm
QN 5070	Size 3-0, 76cm, 1/2 circle RB 25mm
QN 5028	Size 3-0, 76cm, 3/8 circle RC 26mm
QN 5610	Size 2-0, 76cm, 1/2 circle TC 17mm
QN 5670	Size 2-0, 76cm, 1/2 circle TC 25mm
QN 5333	Size 2-0, 76cm, 1/2 circle RB 30mm
QN 5003	Size 2-0, 75cm, 1/2 circle RC 38mm
QN 5036	Size 2-0, 76cm, 3/8 circle RC 45mm
QN 5052	Size 2-0, 76cm, 1/2 circle RC 50mm
QN 5290	Size 2-0, 76cm, Straight cutting 60mm
QN 5017	Size 0, 76cm, 3/8 circle RB 30mm
QN 5334	Size 0, 76cm, 1/2 circle RB 30mm
QN 5004	Size 0, 75cm, 1/2 circle RC 38mm
QN 5037	Size 0, 76cm, 3/8 circle RC 45mm
QN 5005	Size 1, 75cm, 1/2 circle RC 38mm
QN 5062	Size 1, 76cm, 3/8 circle cutting 60mm
QN 5061	Size 4, 76cm, 1/2 circle RB 35mm
QN 5334VS	Size 0, 90cm, 1/2 Circle Round Body
QN 5087A	Size 3-0, 45cm, 1/2 Circle Round Body
QN 5086N	Size 4-0, 76cm, 1/2 Circle Round Body
QN 5070N	Size 3-0, 76cm, 1/2 Circle Round Body
QN 5063	Size 1, 76cm, 1/2 Circle Round Body Heavy
QN 5062N	Size 1, 76cm, 3/8 Circle Cutting 60mm
QN 5002N	Size 3-0, 70cm, 3/8 Circle Cutting
QN 5028A	Size 3-0, 45cm, 3/8 Circle Reverse Cutting
QS 213	Size 2-0, 2x75 cm, without needle
QS 214	Size 0, 2x75 cm, without needle
QS 215	Size 1, 2x75 cm, without needle
QS 212	Size 3-0, 2x75 cm, without needle

Sterility: EO

Packaging: Superior PP box packaging (12 units)

Usage: Single-use only