

## TECHNICAL DATA SHEET



Quickfix™ Skin Stapler  
Disposable surgical tissue closure device

### Technical Specifications

PIN DESCRIPTION: - Length, width and diameter of the pins are determined as follows

Product Code	No. of Pins (+/- 1 Pin)	Pin Shape	Diameter of the pin (mm)	Width of the pin (mm)	Length of the pin (mm)
35W	35	Should be rectangle in shape	0.60 +/- 0.1 mm	4.3 +/- 0.3 mm	7.3 +/- 0.3 mm
45W	45				
55W	55				

**Clinical Benefits:** This easy to use sterile Skin Staplers are used in a variety of surgical procedures for effectively closure of surgical incisions and wounds while minimizing the exposure and operating time and accelerating the rate of closure.

**Sterility:** EO

**Packaging:** Box of 10 sterile units

**Usage:** Single-use only