



TECHNICAL DATA SHEET



SureTack™
Absorbable Tacker
Mesh fixation device



Technical Specifications

Product Code/ Sizes	Shaft Length (mm)	Shaft Dia (mm)	No. of Tacks	Tack Dia (mm)	Tack Height (mm)
QAT 20LN	357.5	5 +/- 0.5	20	5 +/- 0.1	5 +/- 0.1
QAT 30L			30		
Sterility: EO		Packaging: Box of 1 sterile unit		Usage: Single-use only	