



# TECHNICAL DATA SHEET



TiTack™  
Non-absorbable Tacker  
Mesh fixation device



## Technical Specifications

Product Code/ Sizes	Shaft Length (mm)	Shaft Dia (mm)	No. of Tacks	Tack Dia (mm)	Tack Height (mm)
QAT 20LN	357.5	5 +/- 0.5	20	4 +/- 0.1	4 +/- 0.1
QAT 30L			30		
<b>Sterility: EO</b>		<b>Packaging: Box of 1 sterile unit</b>		<b>Usage: Single-use only</b>	